


# HOW IT WORKS

## Battelle CCDS Critical Care Decontamination System™

# BATTELLE




Collect PPE




Package PPE

Items are shipped to the Battelle CCDS™



Items are barcoded to ensure Chain of Custody



HOSPITAL A

HOSPITAL B


PPE is loaded into system to undergo decontamination cycle

Vapor Phase Hydrogen Peroxide Process

TENS OF THOUSANDS of PPE per cycle

Can be processed multiple times without degrading filter performance

PPE is verified to ensure it is FREE OF DECONTAMINANT



PPE is packaged and staged for return

+ Decontaminated PPE is deployed for reuse



Battelle CCDS™ protects healthcare providers, so they can protect us

